# SNAP<sup>®</sup> Residue Tests

Accurate residue detection in 10 minutes or less

## SNAP Test Drug Detection

Test	Drug	Name brand
● <b>New SNAP<sup>®</sup> Beta-Lactam</b>	penicillin G	Albacillin <sup>®</sup> , ALBADRY PLUS <sup>®</sup> , Aqua-Mast II, Benza-Pen, Combicillin <sup>®</sup> , Combi-Pen-48 <sup>®</sup> , Crysticillin, Flo-cillin <sup>®</sup> /Dura-biotic, Hanford/US Vet Han-Pen <sup>®</sup> B (or G), Hanford/US Vet Ultrapen, Hanford/US Vet Masti-Clear <sup>®</sup> , Microcillin, Norocillin <sup>®</sup> , Pen BP-48 <sup>®</sup> , PFI-Pen G, PEN-G MAX <sup>®</sup> , Pro-Pen-G <sup>®</sup> , Procaine G penicillin, Special Formula 17900-Forte <sup>®</sup> Suspension, Quartermaster <sup>®</sup> , Heticin-K
	amoxicillin	Amoxi-Bol, Amoxi-Inject <sup>®</sup> , Amoxi-Mast <sup>®</sup>
	ampicillin	Polyflex <sup>®</sup> , Princillin bolus
	ceftiofur	Excede <sup>®</sup> , Excenel <sup>®</sup> , NAXCEL <sup>®</sup> , Spectramast <sup>®</sup> DC, Spectramast <sup>®</sup> LC
	cephapirin	Cefa-Dri <sup>®</sup> , Cefa-Lak <sup>®</sup> , ToDAY <sup>®</sup> , ToMORROW <sup>®</sup>
● <b>SNAP<sup>®</sup> Tetracycline</b>	tetracycline	Panmycin <sup>®</sup> bolus, Polyotic <sup>®</sup> tablets, Polyotic <sup>®</sup> powder, Tet-Sol <sup>®</sup> , Tetra-Bac <sup>®</sup> , TetraMed <sup>®</sup> , Solu-Tet <sup>®</sup> 324, Tetracycline HCL soluble powder
	chlortetracycline	Aureomycin <sup>®</sup> , Aureo S 700 <sup>®</sup> , CLTC <sup>®</sup> 100 MR, Chlormax <sup>®</sup> , Pennchlor <sup>®</sup>
	oxytetracycline	Agricyl, Agrymycin, Bio-mycin <sup>®</sup> , Geomycin, Liquamycin <sup>®</sup> , Maxim-200 <sup>®</sup> , Medamycin <sup>®</sup> , Neo-Oxy, Neo-Terramycin <sup>®</sup> , Oxy calf bolus, OXTC <sup>®</sup> , Oxy-Tet <sup>®</sup> , Oxytetracycline HCL soluble powder, Oxyshot <sup>®</sup> , Pennox <sup>®</sup> , Promycin <sup>®</sup> , Terramycin <sup>®</sup> , Tetradure <sup>®</sup> , Tetroxy <sup>®</sup> Status SQ, LA-200 <sup>®</sup>
● <b>SNAP<sup>®</sup> Aflatoxin M1</b>	aflatoxin M1	May be found in milk or milk products obtained from livestock that have ingested contaminated feed.
● <b>SNAP<sup>®</sup> Gentamicin</b>	gentamicin	Residues may be found in milk after the intramammary treatment of lactating cows for mastitis.
● <b>SNAP<sup>®</sup> Sulfamethazine</b>	sulfamethazine	Aureo <sup>®</sup> S 700 G, SulfaSure <sup>®</sup> , Sulmet <sup>®</sup> , Sustain III <sup>®</sup>

## SNAP<sup>®</sup> Residue Tests Test Procedure

### Before you begin

- Work space must be clean and free of drug residues.
- Preheat heater block to 45°C ± 5°C (113°F ± 9°F).
- Place SNAP device on heater block; leave on heater block for duration of test.
- Make sure conjugate pellet is at the bottom of the sample tube. If not, tap tube gently to return pellet to the bottom.
- Agitate the milk sample thoroughly.

### Note:

- Accurate incubation times and temperatures are essential.
- Do not use a SNAP test after its expiration date.
- Refer to package insert for complete instructions.

### 1 Prepare sample

Using pipette, draw up milk to indicator line (450 µL ± 50 µL).<sup>†</sup> Add milk to sample tube.

Gently shake tube side-to-side 3–4 times to dissolve reagent pellet.

Place tube on heater block and incubate:

Beta-Lactam	5
Tetracycline	minutes
Aflatoxin M1	2
Gentamicin	minutes
Sulfamethazine	

### 2 Run test

Pour entire contents of sample tube into sample well of SNAP<sup>®</sup> device.

As the edge of the activation circle begins to disappear, press down firmly until you hear a distinct “snap.” If activation takes more than 60 seconds, result may be inaccurate.<sup>‡</sup>

Incubate device:

Beta-Lactam	4
Tetracycline	minutes
Aflatoxin M1	7
Gentamicin	minutes
Sulfamethazine	


### 3 Read results

Read results within 30 seconds after incubation, either visually or using the SNAPshot<sup>®</sup> DSR Reader.

**Negative**  
Sample spot is darker than or equal to control spot. No drug residue detected.

**Positive**  
Sample spot is lighter than control spot. Drug residue detected.

IDEXX SNAP tests are backed by a dedicated team of technical service professionals.

If you have questions, call IDEXX Technical Services at 1-800-321-0207.

<sup>†</sup>NCIMS-certified laboratories must use a calibrated pipette.

<sup>‡</sup>A slow flow rate may be caused by a thick sample or by using less than 450 µL of sample.

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