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Types de globules rouges

- + Sphérocytes (chiens uniquement)
- + Réticulocytes
- + Agglutination (chiens uniquement)
- + Acanthocytes (chiens uniquement)

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Diagnostic Considerations

Acanthocytes indicate alteration of the lipid composition of the RBC membrane or RBC fragmentation and are associated with such disease states as liver/splenic/kidney disease, iron deficiency, disseminated intravascular coagulopathy, cardiac disease and neoplasia.

This platelet estimate incorporates enumeration of individual platelets and platelets within clumps. Platelet count is above 150 K/μL.

Images

RBC/WBC Field of View

Acanthocytes

RBC/WBC (Zoomed in)

a: Results imported from ProCyte.
 b: WBC results imported from ProCyte, the white blood cell differential has been updated based on cellular analysis.

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
Pathology 3/25/26 10:11 AM

Site: Left Ear

Bacteria, Cocci	3-4+ Numerous coccoid-shaped bacteria seen.
Bacteria, Rods	2+ Moderate rod-shaped bacteria seen.
Yeast	2+ Moderate yeast seen.
WBC	Absent
Mites	Absent

Diagnostic Considerations Results: Both bacteria and yeast seen. Consideration: The co-presence of yeast and bacteria should be interpreted in light of clinical signs and degree of burden. Low numbers of bacteria and yeast may represent commensal flora in unaffected ears, early dysbiosis in mildly affected ears, or true infection in clinically affected ears. As the bacteria and yeast burden increases, and if white blood cells are also present, true infection becomes more likely. Next Step(s): Address any primary and predisposing factors—such as atopic dermatitis, endocrinopathies, excessive hair in the ear canals, stenotic ear conformation, increased cerumen production, otic masses, or frequent ear cleaning. In cases of persistent or recurrent infections, especially those with pus or discharge, evaluate the patient for the presence of biofilms, which can make organisms resistant to antimicrobials and require thorough ear cleaning as part of treatment. Use clinical signs, history, and diagnostics to evaluate for deeper involvement of the middle or inner ear. Administer appropriate antimicrobial and anti-inflammatory therapies if warranted based on clinical assessment.

Images

Visit type: **First time**. Case duration: **Chronic**. Clinical signs: Patient presented signs of otitis including:

- Left ear: itchy, pain, erythema (redness), discharge (waxy, watery).

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Pathology 3/26/26 3:16 PM

Site:	Left Ear
Bacteria, Cocci	3-4+ Numerous coccoid-shaped bacteria seen.
Bacteria, Rods	2+ Moderate rod-shaped bacteria seen.
Yeast	1+ Few yeast seen.
WBC	Present
Mites	Absent
Diagnostic Considerations	<p>Results: Both bacteria and yeast seen. Consideration: The co-presence of yeast and bacteria should be interpreted in light of clinical signs and degree of burden. Low numbers of bacteria and yeast may represent commensal flora in unaffected ears, early dysbiosis in mildly affected ears, or true infection in clinically affected ears. As the bacteria and yeast burden increases, and if white blood cells are also present, true infection becomes more likely. Next Step(s): Address any primary and predisposing factors—such as atopic dermatitis, endocrinopathies, excessive hair in the ear canals, stenotic ear conformation, increased cerumen production, otic masses, or frequent ear cleaning. In cases of persistent or recurrent infections, especially those with pus or discharge, evaluate the patient for the presence of biofilms, which can make organisms resistant to antimicrobials and require thorough ear cleaning as part of treatment. Use clinical signs, history, and diagnostics to evaluate for deeper involvement of the middle or inner ear. Administer appropriate antimicrobial and anti-inflammatory therapies if warranted based on clinical assessment.</p> <p>Results: White blood cells present. Consideration: Consider underlying causes of otitis externa—such as atopic dermatitis (food or environmentally triggered), tumor, otitis media, foreign body presence, or contact otitis as from otic cleaners/medications or aggressive mechanical cleaning. Next Step(s): In ears undergoing treatment, persistent inflammation indicates the need to investigate for an underlying cause. Typically, these patients require more intensive/longer duration of treatment and more intensive diagnostics (otic irrigation, advanced imaging, and sometimes otic culture). Address any primary, predisposing, and perpetuating factors as part of the therapeutic regimen.</p>

Images



Visit type: **First time**. Case duration: **Acute**. Clinical signs: Patient presented signs of otitis including:

- Left ear: Itchy, pain, odor, erythema (redness).

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
the presence of biofilms, which can make bacteria resistant to antibiotics and require thorough ear cleaning as part of treatment. Use clinical signs, history, and diagnostics to evaluate for deeper involvement of the middle or inner ear. Administer appropriate antimicrobial and anti-inflammatory therapies based on clinical assessment.

User has indicated that there is purulent discharge in this ear. Consider suppurative otitis: clinically characterized by erythema, ulceration, and a purulent discharge often with a biofilm. These cases usually involve a *Pseudomonas* spp. infection but can rarely be associated with *Staphylococcus* or *Malassezia*. Address potential primary, predisposing, and perpetuating factors such as atopic dermatitis (food or environmentally triggered), tumor, otitis media, foreign body presence, infection and its potential extension to the middle ear, and address biofilm as part of your elected treatment as these protect bacterial colonies from antimicrobial therapy.

Images

Bacteria Assessment Yeast and WBC Assessment (Composite) Yeast and WBC Assessment (Brightfield)

Visit type: First time. Case duration: Chronic. Clinical signs: Patient presented signs of otitis including:
Left ear: pain, itchy, discharge (purulent (pus), bloody)


Pathology 5/14/25 8:02 AM

Source: Right ear

Bacteria, Rods	0	None
Bacteria, Cocci	0	None to trace seen
Yeast	0	None to trace seen
WBC	Absent	
Mites	Present	

Diagnostic Considerations: Results: Otodectes otitis.
Consideration: Any co-presence of bacteria, yeast, and/or white blood cells is likely secondary to ear mite infestation.
Next Step(s): Use an effective acaricide and manage any secondary infections or inflammation.

Images





Bacteria Assessment Yeast and WBC Assessment (Composite) Yeast and WBC Assessment (Brightfield)

Visit type: First time. Case duration: Chronic. Clinical signs: Patient presented signs of otitis including:
Left ear: pain, itchy, discharge (purulent (pus), bloody)

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Source	Right Ear	
Bacteria		
Cocci	3-4+	Numerous coccoid-shaped bacteria seen.
Rods	0	None.
Yeast	1+	Few yeast seen.
WBC	Absent	
Mites	Present	
Results: Otodectes otitis. Consideration: Any co-presence of bacteria, yeast, and/or white blood cells is likely secondary to ear mite infestation. Next Step(s): Use an effective acaricide and manage any secondary infections or inflammation.		
Diagnostic Considerations		
Images		
	Bacteria Assessment	Yeast and WBC Assessment (Composite)
		Yeast and WBC Assessment (Brightfield)

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880012398602
Feline

Bacteria Assessment Yeast and WBC Assessment (Composite) Yeast and WBC Assessment (Brightfield)

Mite Assessment

Page 1 of 1
Tap to enlarge 4 Images in record

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Mites Absent

Diagnostic Considerations

Results: Both rods and cocci observed
Consideration: The co-presence of rods and cocci support bacterial otitis.
Next Step(s): In cases of persistent or recurrent infections, especially those with pus or discharge, evaluate the patient for the presence of biofilms, which can make bacteria resistant to antibiotics and require thorough ear cleaning as part of treatment. Use clinical signs, history, and diagnostics to evaluate for deeper involvement of the middle or inner ear. Administer appropriate antimicrobial and anti-inflammatory therapies based on clinical assessment.

Images

Bacteria Assessment Bacteria Assessment Yeast and WBC Assessment (Brightfield) +7

Source Right Ear

Buxton Brooks
Canine | Poodle | Female | 4 y

Bacteria Assessment Bacteria Assessment Yeast and WBC Assessment (Brightfield)

Yeast and WBC Assessment (Brightfield) Yeast and WBC Assessment (Brightfield) Yeast and WBC Assessment (Brightfield)

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© 2026 IDEXX Laboratories, Inc. Tous droits réservés. • 09-2692856-00
1. Données IDEXX Laboratories, Inc. Westbrook, Maine, États-Unis.