18

Colilert*^{*}-18

06-02027-25

IDEXX

IDX 33/02 - 06/12
IDX 33/01 - 11/09
WATER ANALYSIS METHODS
<http://nf-validation.afnor.org>

The method Colilert*^{*}-18/Quanti-Tray* or Quanti-Tray*/2000 for water analysis is granted NF Validation by AFNOR Certification as an alternative method to the standard ISO 9308-3 for enumeration of *Escherichia coli* β -glucuronidase positive in bathing water, under the Certificate number: IDX 33/02-06/12.

The method Colilert*^{*}-18/Quanti-Tray* for water analysis is granted NF Validation by AFNOR Certification as an alternative method to the standard ISO 9308-1 for enumeration of *Escherichia coli* β -glucuronidase positive and coliform bacteria β -galactosidase positive in drinking water, under the Certificate number: IDX 33/01-11/09.

For more information about end of validity, please refer to the certificate NF Validation available on website mentioned above.

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IDEXX Water Quality Control Laboratory is accredited to ISO/IEC 17025:2017

IDEXX

IDEXX Laboratories, Inc., One IDEXX Drive, Westbrook, Maine 04092 USA
idexx.com/water

For Technical Support, please call:

North/South America: 1 207 556 4496/1 800 321 0207

Europe: 00800 4339 9111

UK: +44 (0) 1638 676800

China: +86 21 61279528

Japan: 03 5301 6800

Australia: 1300 443 399

Colilert® -18 Test Kit

Introduction

Colilert® -18 either simultaneously detects total coliforms and *E. coli*; or fecal coliforms in water. It is based on IDEXX's proprietary Defined Substrate Technology®. When total or fecal coliforms metabolize Colilert-18's DST® nutrient-indicator, ONPG, the sample turns yellow. When *E. coli* metabolize Colilert-18's DST® nutrient-indicator, MUG, the sample also fluoresces. Colilert-18 can simultaneously detect these bacteria at 1 cfu/100 mL within 18 hours even with as many as 2 million heterotrophic bacteria per 100 mL present.

Storage

Store at 2–25°C away from light.

Presence/Absence (P/A) Procedure

1. Add contents of one pack to a 100 mL sample in a sterile, transparent, nonfluorescing vessel.
2. Cap vessel and shake.
3. If sample is not already at 33–38°C, then place vessel in a 35°C waterbath for 20 minutes or, alternatively, a 44.5°C waterbath for 7–10 minutes.
4. Incubate at 35±0.5°C for the remainder of the 18 hours.
5. Read results according to Result Interpretation table below.

Quanti-Tray® Enumeration Procedure

1. Add contents of one pack to a 100 mL room temperature water sample in a sterile vessel.
2. Cap vessel and shake until dissolved.
3. Pour sample/reagent mixture into a Quanti-Tray® or Quanti-Tray®/2000 and seal in an IDEXX Quanti-Tray® Sealer.
4. Place the sealed tray in a 35±0.5°C* (or 44.5±0.2°C for fecal coliforms) incubator for 18 hours (prewarming to 35°C is not required). For incubation in a water bath, submerge the Quanti-Tray, as is, below the water level using a weighted ring.
5. Read results according to the Result Interpretation table below. Count the number of positive wells and refer to the MPN table provided with the trays to obtain a Most Probable Number.

Result Interpretation

Appearance	Result
Less yellow than the comparator ¹ when incubated at 35±0.5°C or 44.5±0.2°C	Negative for total coliforms and <i>E. coli</i> ; Negative for fecal coliforms
Yellow equal to or greater than the comparator when incubated at 35±0.5°C	Positive for total coliforms
Yellow equal to or greater than the comparator when incubated at 44.5±0.2°C	Positive for fecal coliforms
Yellow and fluorescence equal to or greater than the comparator when incubated at 35±0.5°C	Positive for <i>E. coli</i>

- Look for fluorescence with a 6-watt, 365-nm UV light within 5 inches of the sample in a dark environment. Face light away from your eyes and towards the sample.
- Colilert-18 results are to be read after 18 hours of incubation.
- However, if the results are ambiguous to the analyst based on the initial reading, incubate up to an additional four hours (but not to exceed 22 hours total) to allow the color and/or fluorescence to intensify.
- Positives observed before 18 hours and negatives observed after 22 hours are also valid.
- In addition, laboratories may incubate samples for additional time (up to 22 hours total) for their convenience.

Procedural Notes

- A slight tinge may be observed when Colilert-18 is added to the sample.
- If excess foam causes problems while using Quanti-Tray, you may choose to use IDEXX Antifoam Solution (Catalog# WAFDB) or IDEXX 120 ml vessels with Antifoam (Catalog# WW120SBAF-200).
- This insert may not reflect your local regulations. For compliance testing, be sure to follow appropriate regulatory procedures. For example, samples run in other countries should be incubated at 36±2°C for 18–22 hours.
- Colilert-18 can be run in any multiple tube format. *Standard Methods for the Examination of Water and Wastewater*² MPN tables should be used to find Most Probable Numbers (MPNs).
- If a water sample has some background color, compare inoculated Colilert-18 sample to a control blank of the same water sample.
- Colilert-18 can be used for *E. coli* detection (but not coliforms) in marine water. Samples must be diluted at least tenfold. Multiply the MPN value by the dilution factor to obtain the proper quantitative result.
- Use only sterile, nonbuffered, oxidant-free water for dilutions.
- Colilert-18 is a primary water test. Colilert-18 performance characteristics do not apply to samples altered by any pre-enrichment or concentration.
- In samples with excessive chlorine, a blue flash may be seen when adding Colilert-18. If this is seen, consider sample invalid and discontinue testing.
- Aseptic technique should always be followed when using Colilert-18. Dispose of in accordance with Good Laboratory Practices.

Quality Control Procedures—Total Coliform and *E. coli*

1. One of the following quality control procedures is recommended for each lot of Colilert-18:
 - A. IDEXX-QC Coliform and *E. coli*: *Escherichia coli*, *Klebsiella variicola*⁴, and *Pseudomonas aeruginosa*.
 - B. Quanti-Cult[™] *Escherichia coli*, *Klebsiella pneumoniae* and *Pseudomonas aeruginosa*.
 - C. Fill three sterile vessels with 100 mL of sterile non-buffered oxidant-free water and inoculate with a sterile loop of ATCC⁵ strains, *Escherichia coli* ATCC 25922/¹ WDCM 00013 or ATCC 11775/WDCM 00090, *Klebsiella variicola*¹ ATCC 31488/WDCM 00206 or *Klebsiella aerogenes* ATCC 13048/WDCM 00175, and *Pseudomonas aeruginosa* ATCC 10145/WDCM 00024 or ATCC 27853.
 2. Follow the P/A Procedure or Quanti-Tray Enumeration Procedure above.
 3. Results should match the Result Interpretation table above.
- NOTE:** IDEXX internal quality control testing is performed in accordance with ISO 11133:2014. Quality Control Certificates are available at idexx.com/water.

Quality Control Procedures—Fecal Coliform

1. One of the following quality control procedures is recommended for each lot of Colilert-18:
 - A. IDEXX-QC Fecal Coliform⁶: *Escherichia coli* and *Pseudomonas aeruginosa*.
 - B. Quanti-Cult[™] *Escherichia coli* (fecal coliform), *Klebsiella pneumoniae* (fecal coliform) and *Pseudomonas aeruginosa* (non-fecal).
 - C. Fill two sterile vessels with 100 mL of sterile non-buffered oxidant-free water and inoculate with a sterile loop of ATCC⁵ strains, *Escherichia coli* ATCC 11775 (fecal coliform) and *Pseudomonas aeruginosa* ATCC 10145 or 27853 (non-fecal coliform).
2. Follow the Quanti-Tray Enumeration Procedure above.
3. Results should match the Result Interpretation table above.

1. IDEXX P/A Comparator, catalog # WP104; Quanti-Tray Comparator # WQTC, or Quanti-Tray/2000 Comparator # WQ2KC
2. Eaton, AD, Clesceri, LS, Greenberg, AE, Rice, EN. *Standard Methods for the Examination of Water and Wastewater*. American Public Health Association, 2005. Washington, DC.
3. IDEXX QC Coliform and *E. coli*, IDEXX Catalog # UN3373-WQC-TCEC
4. Quanti-Cult cultures—IDEXX catalog # WQIT-1001
5. American Type Culture Collection 1-800-638-6597 atcc.org
6. IDEXX-QC Fecal Coliform—IDEXX Catalog # UN3373-WDC-FC
†. When following AFNOR validated method for drinking water or bathing water testing, place the sealed tray in a 36±2°C incubator for 18 hours (prewarming to 36°C is not required).
‡. *Klebsiella pneumoniae* (ATCC 31488/WDCM 00206) has been renamed to *Klebsiella variicola*.

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© 2019 IDEXX Laboratories, Inc. Alle Rechte vorbehalten. • Patentinformationen: idexx.com/patents.
© 2019 IDEXX Laboratories, Inc. Wszelkie prawa zastrzeżone. • Informacje o patentach: idexx.com/patents.